

Orthopedic




KM39534 : VANGUARD SPLINTER PICK

5" (12.5cm)




KM39536 : POINT ONLY FOR VANGUARD




KM39692 : RULER, INCH & MILLIMETER GRADUATED

6" (152mm), flexible




KM39694 : RULER, INCH & MILLIMETER GRADUATED

8" (203mm), flexible



KM39698 : RULER, INCH & MILLIMETER GRADUATED

12" (305mm), flexible



KM39700 : RULER, INCH & MILLIMETER GRADUATED

20" (508mm), flexible

NO IMAGE
AVAILABLE

KM42618 : PROBE WITH EYE

10" (25.5cm)

NO IMAGE
AVAILABLE

KM42823 : FRAZIER SUCTION TUBE

6 French tip

NO IMAGE
AVAILABLE

KM42824 : FRAZIER SUCTION TUBE

7 French tip

NO IMAGE
AVAILABLE

KM42825 : FRAZIER SUCTION TUBE

8 French tip

NO IMAGE
AVAILABLE

KM42825P : STYLETS ONLY FOR KM 42825 FRAZIER TUBE 8FR

NO IMAGE
AVAILABLE

KM42826 : FRAZIER SUCTION TUBE

10 French tip

NO IMAGE
AVAILABLE

KM42830 : FRAZIER SUCTION TUBE

12 French tip



KM42863 : FRAZIER SUCTION TUBE

Right Angle, 7 French tip, 4 1/2"



KM42864 : FRAZIER SUCTION TUBE

Right Angle, 9 French tip, 4 1/2"



KM42866 : FRAZIER SUCTION TUBE

Right Angle, 11 French tip, 4 1/2"



KM42873 : FRAZIER SUCTION TUBE

Right Angle, 7 French tip, 7 1/2"



KM42874 : FRAZIER SUCTION TUBE

Right Angle, 9 French tip, 7 1/2"



KM42876 : FRAZIER SUCTION TUBE

11 French tip, Working length 7 1/2" (19cm)



KM44030 : APPROXIMATION CLAMP

Straight, 1/2" tip length, 1/2" (1.5cm)

AVAILABLE

NO IMAGE
AVAILABLE

KM44031 : MICRO CLAMP CASE

2 3/4" (70mm) x 1 3/4" (44mm) x 3/8" (9mm)

NO IMAGE
AVAILABLE

KM45208 : SCREW HOLDING SELF-RETAINING FORCEPS 1.5 TO 2.7MM

Angled, 3-3/8" (9 cm)

NO IMAGE
AVAILABLE

KM45210 : SCREW HOLDING SELF-RETAINING FORCEPS 3.5 TO 6.5MM

Angled, 3-3/8" (9 cm)

NO IMAGE
AVAILABLE

KM45302 : BONE HOLDING FORCEPS W/.035 K-WIRE GUIDE

Pointed tips, 6-3/4" (17 cm)

NO IMAGE
AVAILABLE

KM45305 : BONE HOLDING FORCEPS W/.035 K-WIRE GUIDE

Pointed tips, 5" (13 cm)