

## Orthopedic



NO IMAGE  
AVAILABLE

**KM46329 : LAMBOTTE OSTEOTOME**

Curved, 1/2" tip width, 9" (23cm)



NO IMAGE  
AVAILABLE

**KM46330 : LAMBOTTE OSTEOTOME**

Curved, 3/4" tip width, 9" (23cm)



NO IMAGE  
AVAILABLE

**KM46331 : LAMBOTTE OSTEOTOME**

Curved, 1" tip width, 9" (23cm)



NO IMAGE  
AVAILABLE

**KM46332 : LAMBOTTE OSTEOTOME**

Curved, 1 1/4" tip width, 9" (23cm)



NO IMAGE  
AVAILABLE

**KM46333 : LAMBOTTE OSTEOTOME**

Curved, 1 1/2" tip width, 9" (23cm)



NO IMAGE  
AVAILABLE

**KM46335 : LAMBOTTE OSTEOTOME**

Straight, 20mm tip width, 5" (12.5cm)

NO IMAGE  
AVAILABLE

**KM46336 : LAMBOTTE OSTEOTOME**

Straight, 4mm tip width, 5" (12.5cm)

NO IMAGE  
AVAILABLE

**KM46338 : LAMBOTTE OSTEOTOME**

Straight, 6mm tip width, 5" (12.5cm)

NO IMAGE  
AVAILABLE

**KM46340 : LAMBOTTE OSTEOTOME**

Straight, 8mm tip width, 5" (12.5cm)

NO IMAGE  
AVAILABLE

**KM46342 : LAMBOTTE OSTEOTOME**

Straight, 10mm tip width, 5" (12.5cm)

NO IMAGE  
AVAILABLE

**KM46343 : LAMBOTTE OSTEOTOME**

Straight, 14mm tip width, 5" (12.5cm)

NO IMAGE  
AVAILABLE

**KM46344 : LAMBOTTE OSTEOTOME**

Straight, 12mm tip width, 5" (12.5cm)

NO IMAGE  
AVAILABLE

**KM46345 : LAMBOTTE OSTEOTOME**

Straight, 15mm tip width, 5" (12.5cm)



**KM46346 : LAMBOTTE OSTEOTOME**

Straight, 18mm tip width, 5" (12.5cm)



**KM46347 : LAMBOTTE OSTEOTOME**

Straight, 19mm tip width, 5" (12.5cm)



**KM46348 : LAMBOTTE OSTEOTOME**

Straight, 4mm tip width, with calibration lines, 7" (18cm)



**KM46349 : LAMBOTTE OSTEOTOME**

Straight, 5mm tip width, with calibration lines, 7" (18cm)



**KM46350 : LAMBOTTE OSTEOTOME**

Straight, 6mm tip width, with calibration lines, 7" (18cm)



**KM46351 : LAMBOTTE OSTEOTOME**

Straight, 10mm tip width, with calibration lines, 7" (18cm)



**KM46352 : LAMBOTTE OSTEOTOME**

Straight, 12mm tip width, with calibration lines, 7" (18cm)

AVAILABLE

NO IMAGE  
AVAILABLE

**KM46353 : LAMBOTTE OSTEOTOME**

Straight, 16mm tip width, with calibration lines, 7" (18cm)

NO IMAGE  
AVAILABLE

**KM46354 : LAMBOTTE OSTEOTOME**

Straight, 19mm tip width, with calibration lines, 7" (18cm)

NO IMAGE  
AVAILABLE

**KM46355 : LAMBOTTE OSTEOTOME**

Straight, 22mm tip width, with calibration lines, 7" (18cm)

NO IMAGE  
AVAILABLE

**KM46356 : LAMBOTTE OSTEOTOME**

Straight, 25mm tip width, with calibration lines, 7" (18cm)

NO IMAGE  
AVAILABLE

**KM46358 : LAMBOTTE OSTEOTOME**

Straight, 8mm tip width, with calibration lines, 7" (18cm)