

Orthopedic




KM46459 : SMITH PETERSON OSTEOTOME

Curved, 1/2" tip width, 8" (20.5cm)




KM46460 : SMITH PETERSON OSTEOTOME

Curved, 5/8" tip width, 8" (20.5cm)




KM46461 : SMITH PETERSON OSTEOTOME

Curved, 3/4" tip width, 8" (20.5cm)




KM46462 : SMITH PETERSON OSTEOTOME

Curved, 1" tip width, 8" (20.5cm)



KM46463 : SMITH-PETERSON OSTEOTOME

Curved, 1 1/4" tip width, 8" (20.5cm)



KM46466 : SMITH PETERSON GOUGE

Straight, 1/4" distal tip, 8" (20.3cm)

NO IMAGE
AVAILABLE

KM46467 : SMITH PETERSON GOUGE

Straight, 3/8" distal tip, 8" (20.3cm)

NO IMAGE
AVAILABLE

KM46468 : SMITH PETERSON GOUGE

Straight, 1/2" distal tip, 8" (20.3cm)

NO IMAGE
AVAILABLE

KM46476 : SMITH PETERSON GOUGE

Straight, 5/8" distal tip, 8" (20.3cm)

NO IMAGE
AVAILABLE

KM46477 : SMITH PETERSON GOUGE

Straight, 3/4" distal tip, 8" (20.3cm)

NO IMAGE
AVAILABLE

KM46478 : SMITH PETERSON GOUGE

Straight, 1" distal tip, 8" (20.3cm)

NO IMAGE
AVAILABLE

KM46479 : SMITH PETERSON GOUGE

Straight, 1 1/4" distal tip, 8" (20.3cm)

NO IMAGE
AVAILABLE

KM46487 : SMITH PETERSON GOUGE

Curved, 1/4" distal tip, 8" (20.3cm)

NO IMAGE
AVAILABLE

KM46488 : SMITH PETERSON GOUGE

Curved, 3/8" distal tip, 8" (20.3cm)

NO IMAGE
AVAILABLE

KM46489 : SMITH PETERSON GOUGE

Curved, 1/2" distal tip, 8" (20.3cm)

NO IMAGE
AVAILABLE

KM46490 : SMITH PETERSON GOUGE

Curved, 5/8" distal tip, 8" (20.3cm)

NO IMAGE
AVAILABLE

KM46491 : SMITH PETERSON GOUGE

Curved, 3/4" distal tip, 8" (20.3cm)

NO IMAGE
AVAILABLE

KM46493 : SMITH PETERSON GOUGE

Curved, 1" distal tip, 8" (20.3cm)

NO IMAGE
AVAILABLE

KM46495 : SMITH-PETERSON GOUGE

Curved, 1 1/4" distal tip, 8" (20.3cm)

NO IMAGE

KM46497 : HOKE OSTEOTOME

Straight, 1/8" (3mm) tip width, 5 1/4" (13.5cm)

AVAILABLE

NO IMAGE
AVAILABLE

KM46498 : HOKE OSTEOTOME

Straight, 3/16" (4mm) tip width, 5 1/4" (13.5cm)

NO IMAGE
AVAILABLE

KM46499 : HOKE OSTEOTOME

Straight, 1/4" (6mm) tip width, 5 1/4" (13.5cm)

NO IMAGE
AVAILABLE

KM46500 : HOKE OSTEOTOME

Straight, 5 16" (8mm) tip width, 5 1/4" (13.5cm)

NO IMAGE
AVAILABLE

KM46501 : HOKE OSTEOTOME

Straight, 3/8" (10mm) tip width, 5 1/4" (13.5cm)

NO IMAGE
AVAILABLE

KM46502 : HOKE OSTEOTOME

Straight, 1/2" (13mm) tip width, 5 1/4" (13.5cm)