

Orthopedic



NO IMAGE
AVAILABLE

KM47784 : LOVE-GRUENWALD IVD RONGEUR

Bite down, 3 x 10mm, Working length 7" (18cm)



NO IMAGE
AVAILABLE

KM47800 : SPURLING IVD RONGEUR

Straight, 4 x 10mm, Working length 5" (12.5cm)



NO IMAGE
AVAILABLE

KM47802 : SPURLING IVD RONGEUR

Bite up, 4 x 10mm, Working length 5" (12.5cm)



NO IMAGE
AVAILABLE

KM47804 : SPURLING IVD RONGEUR

Bite down, 4 x 10mm, Working length 5" (12.5cm)



NO IMAGE
AVAILABLE

KM47840 : SPURLING IVD RONGEUR

Straight, 4 x 10mm, Working length 7" (18cm)



NO IMAGE
AVAILABLE

KM47842 : SPURLING IVD RONGEUR

Bite up, 4 x 10mm, Working length 7" (18cm)

NO IMAGE
AVAILABLE

KM47844 : SPURLING IVD RONGEUR

Bite down, 4 x 10mm, Working length 7" (18cm)

NO IMAGE
AVAILABLE

KM47864 : FERRIS-SMITH IVD RONGEUR

Ring handle, Straight, 3 x 10mm, Working length 7" (18cm)

NO IMAGE
AVAILABLE

KM47866 : FERRIS-SMITH IVD RONGEUR

Ring handle, Bite up, 3 x 10mm, Working length 7" (18cm)

NO IMAGE
AVAILABLE

KM47868 : FERRIS-SMITH IVD RONGEUR

Ring handle, Bite down, 3 x 10mm, Working length 7" (18cm)

NO IMAGE
AVAILABLE

KM47870 : FERRIS-SMITH IVD RONGEUR

Ring handle, Straight, 4 x 10mm, Working length 7" (18cm)

NO IMAGE
AVAILABLE

KM47872 : FERRIS-SMITH IVD RONGEUR

Ring handle, Bite up, 4 x 10mm, Working length 7" (18cm)

NO IMAGE
AVAILABLE

KM47873 : FERRIS-SMITH IVD RONGEUR

Ring handle, Bite up, 6 x 12mm, Working length 7" (18cm)

NO IMAGE
AVAILABLE

KM47874 : FERRIS-SMITH IVD RONGEUR

Ring handle, Bite down, 4 x 10mm, Working length 7" (18cm)

NO IMAGE
AVAILABLE

KM47875 : FERRIS-SMITH IVD RONGEUR

Ring handle, Bite down, 6 x 12mm, Working length 7" (18cm)

NO IMAGE
AVAILABLE

KM47878 : FERRIS-SMITH IVD RONGEUR

Ring handle, Straight, 6 x 10mm, Working length 7" (18cm)

NO IMAGE
AVAILABLE

KM47879 : SPURLING KERRISON LAMINECTOMY RONGEUR

Angled forward, 1mm, Working length 7" (18cm)

NO IMAGE
AVAILABLE

KM47907 : FIXATION PIN

NO IMAGE
AVAILABLE

KM47908 : Mini Quick Coupling Rounded Phenolic Handle

For holding quick coupling drills, countersinks, taps and screwdrivers. For mini screw set. Features a rounded phenolic handle. 4-3/8" (11 cm). *Note: image shows KM47902 Mini Quick Coupling Handle without rounded phenolic...*

NO IMAGE
AVAILABLE

KM47919 : PLATE BENDING IRON

for 3.5mm and 4.5mm plates

AVAILABLE

NO IMAGE
AVAILABLE

KM47927 : PLATE BENDING PLIERS

For reconstruction plates, 9 3/8" (23.5cm)

NO IMAGE
AVAILABLE

KM47930 : SCREW HOLDING SELF-RETAINING FORCEPS

Curved, 3-1/2" (9 cm)

NO IMAGE
AVAILABLE

KM47932 : PLATE CUTTING FORCEPS UP TO 1MM

Straight, 9-1/4" (24 cm)

NO IMAGE
AVAILABLE

KM47996 : SPURLING KERRISON LAMINECTOMY RONGEUR

Angled forward, 2mm, Working length 7 1/8" (18cm)

NO IMAGE
AVAILABLE

KM47997 : SPURLING KERRISON LAMINECTOMY RONGEUR

Angled forward, 3mm, Working length 7 1/8" (18cm)