

Orthopedic



NO IMAGE
AVAILABLE

KM47998 : SPURLING KERRISON LAMINECTOMY RONGEUR

Angled forward, 4mm, Working length 7 1/8" (18cm)



NO IMAGE
AVAILABLE

KM47999 : SPURLING KERRISON LAMINECTOMY RONGEUR

Angled forward, 5mm, Working length 7 1/8" (18cm)



NO IMAGE
AVAILABLE

KM48010 : FERRIS-SMITH KERRISON LAMINECTOMY RONGEUR

1mm



NO IMAGE
AVAILABLE

KM48030 : FORAMINOTOMY PUNCH

Curved up, 2mm, Working length 8" (20.5cm)



NO IMAGE
AVAILABLE

KM48032 : FORAMINOTOMY PUNCH

Curved up, 4mm, Working length 8" (20.5cm)



NO IMAGE
AVAILABLE

KM48033 : FORAMINOTOMY PUNCH

Curved up, 3mm, Working length 8" (20.5cm)

NO IMAGE
AVAILABLE

KM48041 : FERRIS-SMITH KERRISON RONGEUR

Bite down 90, 6mm, Working length 7" (18cm)

NO IMAGE
AVAILABLE

KM48242 : NEEDLE NOSE PLIERS

Delicate, 5 3/8" (13.5cm)

NO IMAGE
AVAILABLE

KM48243 : FLAT NOSE PLIERS

Combination plier and cutter, 6 1/2" (16.5cm)

NO IMAGE
AVAILABLE

KM48245 : NEEDLE NOSE PLIERS

5 3/8" (13.5cm)

NO IMAGE
AVAILABLE

KM48246 : FLAT NOSE PLIERS

Broad, 4mm tip, 5 3/8" (13.5cm)

NO IMAGE
AVAILABLE

KM48248 : FLAT NOSE PLIERS

Delicate, 1.5mm tip, 6" (15cm)

NO IMAGE
AVAILABLE

KM48250 : LONG JAW NEEDLE NOSE PLIERS

Parallel grooved 2" (50mm) jaws, 1.5mm tip, 7" (18cm)



KM48253 : NARROW NOSE WIRE PULLING FORCEPS

3mm tip, 7 1/2" (19cm)



KM48254 : FLAT NOSE K-WIRE PLIERS

3mm tip, 6 3/4" (17cm)



KM48256 : FLAT NOSE K-WIRE PLIERS

3mm tip, 6 3/4" (17cm)



KM48257 : K-WIRE EXTRACTION PLIERS W/SCREW LOCK

3mm tip, 5" (13cm)



KM48258 : K-WIRE EXTRACTION PLIERS W/SCREW LOCK

4mm tip, 7" (18cm)



KM48259 : LINEMAN PLIERS

Side cutting, maximum capacity .062" (1.6mm), 8 1/2" (21.5cm)



KM48262 : SLIP JOINT PLIERS

8" (20.5cm)

AVAILABLE

NO IMAGE
AVAILABLE

KM48269 : WIRE BENDING PLIERS

5 7/8" (15cm)

NO IMAGE
AVAILABLE

KM48270 : KIRSCHNER BOW SMALL

5 3/4" (14.5cm)

NO IMAGE
AVAILABLE

KM48272 : KIRSCHNER BOW LARGE

7 1/2" (19.0cm)

NO IMAGE
AVAILABLE

KM48274 : KIRSCHNER BOW EXTRA LONG

9 5/8" (24.5cm)

NO IMAGE
AVAILABLE

KM48276 : GUIDEWIRE W/BALL 2 MM x 60 CM