

## Orthopedic



### **KM49362 : COBB OSTEOTOME**

Curved, 1" tip width, 11" (28cm)




### **KM49363 : COBB OSTEOTOME**

Curved, 1 1/4" tip width, 11" (28cm)




### **KM49574 : SCOVILLE CURETTE**

Straight, 10" (25.5 cm)




### **KM49576 : SCOVILLE CURETTE**

Down angle, 10" (25.5 cm)



### **KM49578 : SCOVILLE CURETTE**

Down angle, 10" (25.5 cm)



### **KM49580 : EPSTEIN CURETTE**

Reverse angle, Size 000, 8" (20.5 cm)

NO IMAGE  
AVAILABLE

**KM49582 : EPSTEIN CURETTE**

Reverse angle, Size 00, 8" (20.5 cm)

NO IMAGE  
AVAILABLE

**KM49584 : EPSTEIN CURETTE**

Reverse angle, Size 0, 8" (20.5 cm)

NO IMAGE  
AVAILABLE

**KM49628 : ALLIS WILLAUER FORCEPS**

5 x 6 teeth, 9-7/8" (25 cm)

NO IMAGE  
AVAILABLE

**KM49630 : INTERCHANGEABLE OSTEOTOME, CHISEL AND  
GOUGE SET**

NO IMAGE  
AVAILABLE

**KM49640 : GOUGE BLADE FOR KM49630 (INTERCHANGEABLE  
OSTEOTOME , CHISEL AND GOUGE SET)**

1" (50mm radius) distal tip

NO IMAGE  
AVAILABLE

**KM49644 : GOUGE BLADE FOR KM49630 (INTERCHANGEABLE  
OSTEOTOME , CHISEL AND GOUGE SET)**

1" (60mm radius) distal tip

NO IMAGE  
AVAILABLE

**KM49648 : GOUGE BLADE FOR KM49630 (INTERCHANGEABLE  
OSTEOTOME , CHISEL AND GOUGE SET)**

1" (70mm radius) distal tip




**KM49652 : Chisel Blade For KM49630 (Interchangeable Osteotome , Chisel and Gouge Set)**

Straight, 5mm tip




**KM49656 : Chisel Blade For KM49630 (Interchangeable Osteotome , Chisel and Gouge Set)**

Straight, 10mm tip




**KM49660 : Chisel Blade For KM49630 (Interchangeable Osteotome , Chisel and Gouge Set)**

Straight, 16mm tip




**KM49664 : Chisel Blade For KM49630 (Interchangeable Osteotome , Chisel and Gouge Set)**

Straight, 25mm tip




**KM49668 : OSTEOTOME BLADE FOR KM49630 (INTERCHANGEABLE OSTEOTOME , CHISEL AND GOUGE SET)**

Straight, 5mm tip width, 8cm



**KM49672 : OSTEOTOME BLADE FOR KM49630 (INTERCHANGEABLE OSTEOTOME , CHISEL AND GOUGE SET)**

Straight, 10mm tip width, 8cm



**KM49676 : OSTEOTOME BLADE FOR KM49630 (INTERCHANGEABLE OSTEOTOME , CHISEL AND GOUGE SET)**

Straight, 16mm tip width, 8cm

AVAILABLE

NO IMAGE  
AVAILABLE

**KM49680 : OSTEOTOME BLADE FOR KM49630  
(INTERCHANGEABLE OSTEOTOME , CHISEL AND GOUGE SET)**

Straight, 25mm tip width, 8cm

NO IMAGE  
AVAILABLE

**KM49706 : ROD HOLDER 5.0MM SMOOTH**

8 1/4" (21.0cm)

NO IMAGE  
AVAILABLE

**KM49724 : K-WIRE&PLATE BENDER 3-PRONG MAX.1.1MM**

5" (12.5cm)

NO IMAGE  
AVAILABLE

**KM49734 : K-WIRE SIDE CUTTER**

Welded carbide jaws, max. capacity: .045" (1.1mm), 6" (15cm)

NO IMAGE  
AVAILABLE

**KM49736 : KMEDIC SIDE CUTTER**

Welded carbide jaws, max cap: 3/32" (2.4mm), 9" (23cm)