

Orthopedic

NO IMAGE
AVAILABLE

KM49738 : KMEDIC END CUTTER

Welded carbide jaws, max cap: 3/32" (2.4mm), 8 1/2" (21.5cm)

NO IMAGE
AVAILABLE

KM49740 : KMEDIC WIRE CUTTER WITH SAFETY INSERTS

Max. capacity: .062" , 6 7/8" (17.5cm)

NO IMAGE
AVAILABLE

KM49740P : SILICON INSERTS FOR KM 49-740

NO IMAGE
AVAILABLE

KM49750 : BONE CHIP PACKER

Head diameter 3mm x 16mm, 9 7/8" (25.0cm)

NO IMAGE
AVAILABLE

KM49751 : BONE CHIP PACKER

Head diameter 2.4mm x 9.5mm, 9 7/8" (25.0cm)

NO IMAGE
AVAILABLE

KM49752 : BONE CHIP PACKER

Head diameter 6.4mm x 9.5mm, 9 7/8" (25.0cm)

NO IMAGE
AVAILABLE

KM49753 : BONE CHIP PACKER

Head diameter 6.4mm x 19mm, 9 7/8" (25.0cm)

NO IMAGE
AVAILABLE

KM49754 : BONE CHIP PACKER

Head diameter 1/4" (0.5cm), 9 7/8" (25.0cm)

NO IMAGE
AVAILABLE

KM49760 : COLLIN MALLET

Solid, Head diameter 30.0mm, Weight 7 oz (198 grams); Head weight: 5 oz (142 grams), 8" (20.5cm)

NO IMAGE
AVAILABLE

KM49761 : LUCAE BONE MALLET

Solid, Head diameter 19.0mm, Weight 7 oz (198 grams); Head weight: 3 oz (85 grams), 7 1/2" (19.0cm)

NO IMAGE
AVAILABLE

KM49762 : ORTHOPEDIC MALLET, HEAVYWEIGHT

Solid, Head diameter 38.0mm, Weight 32 oz (907 grams); Head weight: 21 oz (595 grams), 7 1/2" (19.0cm)

NO IMAGE
AVAILABLE

KM49763 : ORTHOPEDIC MALLET, LIGHTWEIGHT

Solid, Head diameter 30.5mm, Weight 18 oz (510 grams); Head weight: 8 oz (227 grams), 7" (18.0cm)

NO IMAGE
AVAILABLE

KM49764 : NARROW TIP MALLET

Head diameter 18 mm x 18mm (square contact surface), Weight 7 oz (198 grams); Head weight: 3 oz (85 grams), 8" (20.5cm)



KM49765 : MINI MALLET

Solid, Head diameter 20.0mm, Weight 4 oz (113 grams); Head weight: 3 oz (85 grams), 7" (18.0cm)



KM49766 : BERGMAN MALLET

Stainless Steel and Aluminum, Head diameter 30.0mm, Weight 15 oz (425 grams); Head weight: 9 oz (255 grams), 8 7/8" (22.5cm)




KM49767 : BERGMAN MALLET

Stainless Steel and Aluminum, Head diameter 45.0mm, Weight 17 oz (482 grams); Head weight: 11 oz (312 grams), 9 1/2" (24.0cm)




KM49780 : IN SITU ROD BENDER 4.76MM RIGHT

13" (33.0cm)




KM49781 : IN SITU ROD BENDER 4.76MM LEFT

13" (33.0cm)



KM49782 : IN SITU ROD BENDER 6.35MM RIGHT

13" (33.0cm)



KM49783 : IN SITU ROD BENDER 6.35MM LEFT

13" (33.0cm)

AVAILABLE

NO IMAGE
AVAILABLE

KM49800 : KMEDIC DOUBLE ACTION PIN CUTTER

Max. capacity" 1/8", 9 1/2" (24cm)

NO IMAGE
AVAILABLE

KM49802 : K-WIRE BENDER 3 PRONG MAX.1.6MM

5 3/4" (14.5cm)

NO IMAGE
AVAILABLE

KM49810 : WIRE TWISTER JET 10MM JAW 11"

11" (28.0cm)

NO IMAGE
AVAILABLE

KM49812 : ROD BENDER FOR 4 & 3MM RODS

7 1/4" (18.5cm)

NO IMAGE
AVAILABLE

KM49832 : ROD HOLDER 6.35MM

7 1/4" (18.5cm)