

Orthopedic



NO IMAGE
AVAILABLE

KM48354 : SCREWDRIVER

10" (25.5cm), Single slot



NO IMAGE
AVAILABLE

KM35460 : Super-Cut Iris Scissors

Curved, 4 1/2" (11.5cm)



NO IMAGE
AVAILABLE

KM44028 : MICRO CLIP APPLICATOR

5 1/2" (14.0cm)



NO IMAGE
AVAILABLE

KM48052 : DRILL BIT QUICK COUPLING

1.5mm dia., 18mm, 3 3/8" (8.5cm), 60mm working length



NO IMAGE
AVAILABLE

KM52158 : RIGID RAKE RETRACTOR SKIN HOOK

Sharp, double-pronged, 6-1/2" (16.5 cm)



NO IMAGE
AVAILABLE

KM46917 : SUTURE PASSER W/CROCHET HOOK

Curved, crochet hook, 9-1/4" (23.5 cm)

NO IMAGE
AVAILABLE

KM48161 : FRAZIER HOOK

Blunt, curved, 2.5 mm diam., single prong, 5" (12.7 cm)

NO IMAGE
AVAILABLE

KM55340 : RHOTON MICRO DISSECTION HOOK

Sharp, 45 shank, 7-1/2" (19 cm)

NO IMAGE
AVAILABLE

KM54330 : GRAEFE STRABISMUS HOOK, SMALL

5-1/8" (13.1 cm)

NO IMAGE
AVAILABLE

KM52164 : RIGID RAKE RETRACTOR SKIN HOOK

Blunt, single prong, 6-1/2" (16.5 cm)

NO IMAGE
AVAILABLE

KM52176 : FLEXIBLE RAKE RETRACTOR SKIN HOOK

Sharp, double-pronged, 6-1/4" (15.9 cm)

NO IMAGE
AVAILABLE

KM52776 : GELPI RETRACTOR

3 1/2" (9cm)

NO IMAGE
AVAILABLE

481250 : Love Nerve Retractor

Blade 4 x 6mm, straight, 8" (20.0cm)

NO IMAGE
AVAILABLE

KM48352 : SCREWDRIVER

3.5mm hexagonal dia., 9 7/8" (25cm)

NO IMAGE
AVAILABLE

KM52156 : RIGID RAKE RETRACTOR SKIN HOOK

Sharp, single prong, 6-1/2" (16.5 cm)

NO IMAGE
AVAILABLE

KM52778 : GELPI RETRACTOR

4 3/8" (11cm)

NO IMAGE
AVAILABLE

KM55339 : RHOTON MICRO DISSECTION HOOK

Blunt, 90 shank, 7-1/2" (19 cm)

NO IMAGE
AVAILABLE

KM48368 : HEXAGONAL SCREWDRIVER

2.5mm with notch dia., 8 1/4" (21cm)

NO IMAGE
AVAILABLE

KM55420 : Blade Retractor

18mm width blade, 7" (18.0cm)

NO IMAGE
AVAILABLE

KM58083 : GILLIES HOOK, LARGE

7" (17.8 cm)

AVAILABLE

NO IMAGE
AVAILABLE

KM48067 : DRILL BIT QUICK COUPLING

4.5mm dia., 50mm, 5 3/4" (14.5cm), 120mm working length

NO IMAGE
AVAILABLE

KM48358 : SCREWDRIVER

6mm dia., 10" (25.5cm)

NO IMAGE
AVAILABLE

KM51080 : Israel Retractor

4-pronged, 1 5/8" (4.0cm) width blade, 9 1/2" (24.0cm)

NO IMAGE
AVAILABLE

KM48210 : Aufranc Cobra Retractor

hollow handle, serrated blunt 10mm tip, 1 1/4" blade width, 10 1/2" (26.5cm)

NO IMAGE
AVAILABLE

KM47555 : KLEINERT-KUTZ BONE CUTTING FORCEPS

Straight, delicate, 5-5/8" (15 cm)