

Orthopedic

NO IMAGE
AVAILABLE**KM55401 : CLOWARD-STYLE MALLET**

Solid, Head diameter 25mm;ball pean: 22mm, Weight 16 oz (454 grams); Head weight: 8 oz (227 grams), 9" (23.0cm)

NO IMAGE
AVAILABLE**KM55498 : D'Errico-Adson Retractor**

4 X 4-pronged, 3/4" (2.0cm) depth blade, 7 3/8" (19.0cm)

NO IMAGE
AVAILABLE**KM55510 : ENDO SCOVILLE CURETTE**

Straight tip, 3-1/2" long, 9 x 15 mm cup, 10" shaft length, 6 1/2" handle length, 20" (51 cm) overall length

NO IMAGE
AVAILABLE**KM55523 : ENDO BRUN CURETTE**

Straight, 15 X 18mm, 17" (43 cm), 9.3", 5"

NO IMAGE
AVAILABLE**KM55533 : ENDO BRUN CURETTE**

Angled, 40, 9 X 15mm, 17" (43 cm), 9.3", 5"

NO IMAGE
AVAILABLE**KM55556 : ENDO BRUN CURETTE**

Angled, 30, 2 X 4mm, 20" (51 cm), 10", 6 1/2"

NO IMAGE
AVAILABLE

KM55566 : ENDO CONE RING CURETTE

Angled, 30, 3 X 6mm, 20" (51 cm), 10", 6 1/2"

NO IMAGE
AVAILABLE

KM55586 : ENDO EPSTEIN CURETTE

Reverse angle, 6 X 9mm, 20" (51 cm), 10", 6 1/2"

NO IMAGE
AVAILABLE

KM55593 : ENDO EPSTEIN CURETTE

Reverse angle, 3 X 6mm, 17" (43 cm), 9.3", 5"

NO IMAGE
AVAILABLE

KM55612 : ENDOSCOPIC COBB ELEVATOR

Sharp, 3/8" (10mm) tip width, 3 1/2" tip length, 10" shaft length, 6 1/2" handle length, .45" (11.5mm shaft dia., 20" (51cm)

NO IMAGE
AVAILABLE

KM55614 : ENDOSCOPIC COBB ELEVATOR

.490" (13mm) tip width, 20" (51cm)

NO IMAGE
AVAILABLE

KM55620 : ENDOSCOPIC COBB ELEVATOR

Sharp, 1/4" (6mm) tip width, 2.7" tip length, 9.3" shaft length, 5" handle length, .285" (7.2mm) shaft dia., 17" (43cm)

NO IMAGE
AVAILABLE


KM55640 : ENDOSCOPIC RIB ELEVATOR

Sharp, Angled, 45, 10mm, 17" (43cm)




KM55662 : NERVE PROBE

Angled, 45, 15" (38cm)




KM55664 : NERVE PROBE

Angled, 90, 15" (38cm)




KM55690 : SCHLESINGER IVD RONGEUR

3mm




KM55691 : SCHLESINGER IVD RONGEUR

Straight, 5mm, 11 3/4" (30cm), Working length 9" (23cm)




KM55763 : FERRIS-SMITH ENDOSCOPIC IVD RONGEUR

Ring handle, Straight, 3.25", 4 x 10mm, 12mm shaft diameter, 20" (51cm), Working length 14" (35.5cm)



KM55777 : FERRIS-SMITH ENDOSCOPIC IVD RONGEUR

Ring handle, Straight, 3.25", 8 x 12mm, 12mm shaft diameter, 20" (51cm), Working length 14" (35.5cm)



KM55782 : FERRIS-SMITH ENDOSCOPIC KERRISON RONGEUR

Ring handle, Bite down 90, 3.25", 4mm, 12mm shaft diameter, 20" (51cm), Working length 14" (35.5cm)

AVAILABLE

NO IMAGE
AVAILABLE

KM55790 : FERRIS-SMITH ENDOSCOPIC KERRISON RONGEUR

Ring handle, Angled forward 40, 3.25", 12mm shaft diameter, 20" (51cm), Working length 14" (35.5cm)

NO IMAGE
AVAILABLE

KM55825 : FERRIS SMITH KERRISON LAMINECTOMY RONGEUR

Ring handle, 5mm

NO IMAGE
AVAILABLE

KM55842 : PEAPOD IVD RONGEUR

Bite up, 2 x 6mm, Working length 5 1/2" (14cm)

NO IMAGE
AVAILABLE

KM55856 : OLDBERG IVD RONGEUR

Straight round, 6mm, Working length 7" (18cm), Round Cup

NO IMAGE
AVAILABLE

KM55858 : DECKER IVD RONGEUR

Straight, 2 x 6mm, Working length 6" (15cm)