

Orthopedic



NO IMAGE
AVAILABLE

KM71895 : STEINMAN PIN

5/32" (4.0mm) Single Diamond, 3 Shank End with suture passer hole, 9", BX/6



NO IMAGE
AVAILABLE

KM71898 : STEINMAN PIN

.177" (4.5mm) Single Diamond, 3 Shank End with suture passer hole, 9", BX/6



NO IMAGE
AVAILABLE

KM98150 : HOKE SMARTRACK OSTEOTOME

- Nylon coated brackets help avoid metal-to-metal contact.
- Plastic plate on underside of top bracket helps protect sharp tips from being damaged.
- Side arm stays locked to help prevent instruments from falling out, even...



NO IMAGE
AVAILABLE

KM98160 : BRUN SMARTRACK FOR 7" CURETTES

Aluminum and Stainless Steel



NO IMAGE
AVAILABLE

KM98170 : COBB SMARTRACK FOR 11" INSTRUMENTS

Aluminum and Stainless Steel



NO IMAGE
AVAILABLE

KM98180 : HIBBS&BRUN SMARTRACK 9" INSTRUMENT

Aluminum and Stainless Steel

NO IMAGE
AVAILABLE

KM98190 : SMARTRACK LAMBOTTE FOR 9" OSTEOTOMES

Aluminum and Stainless Steel

NO IMAGE
AVAILABLE

KM98200 : BONE MILL WITH CRANK & TABLE CLAMP

NO IMAGE
AVAILABLE

KM99100 : PHENOLIC MATERIAL 1000 MM X 32 MM

NO IMAGE
AVAILABLE

OS100199 : TL SET, CASE ASSEMBLY

Case components only. Set consists of: OS100201, OS100202, OS100209 and OS100210.

NO IMAGE
AVAILABLE

P12838 : OCHSNER TROCAR, 12 FR

For 8 Fr catheter

NO IMAGE
AVAILABLE

P14279 : CRILE NERVE HOOK

7 mm blunt tip, 8" (20 cm)

NO IMAGE
AVAILABLE

P23450 : ALLERDYCE ELEVATOR AND DISSECTOR

Double-ended, 8 1/2" (21.5cm)

NO IMAGE
AVAILABLE

P23483 : Smillie Knee Retractor

19 x 30mm, 5 1/2" (14cm)

NO IMAGE
AVAILABLE

P23558 : Kolbel Retractor Complete

w/ 3 pairs of blades

NO IMAGE
AVAILABLE

P23560 : Kolbel Self Retaining Retractor Snap Closure, Frame Only

NO IMAGE
AVAILABLE

P23567 : Kolbel Retractor Blades, Extra Long

pair of 36 x 68mm blades

NO IMAGE
AVAILABLE

P23568 : Kolbel Retractor Blades Set

pair of 36 x 85mm blades

NO IMAGE
AVAILABLE

P23583 : Goldstein Glenoid Neck Retractor

18mm tip, 8 1/2" (21.5cm)

NO IMAGE
AVAILABLE

P23584 : Deltoid Retractor

30mm blade, 8" (20.5cm)

AVAILABLE

NO IMAGE
AVAILABLE

P23586 : Anterior Glenoid Neck Retractor

25mm wide tip, 11" (28cm)

NO IMAGE
AVAILABLE

P23593 : Kolbel Glenoid Retractor

15mm, 11 3/4" (30cm)

NO IMAGE
AVAILABLE

P23594 : Kolbel Glenoid Retractor

2 prong, lightly curved, 15mm, 11 3/4" (30cm)

NO IMAGE
AVAILABLE

P23595 : Kolbel Glenoid Retractor

20mm, 11 3/4" (30cm)

NO IMAGE
AVAILABLE

P23596 : Kolbel Glenoid Retractor

2 prong, lightly curved, 20mm, 11 3/4" (30cm)