

Orthopedic

NO IMAGE
AVAILABLE

P22280 : MCCULLOCH TYPE HOOK BLADE

1 1/2" (4cm)

NO IMAGE
AVAILABLE

P23999 : QUICK COUPLING T-HANDLE, SMALL

3 1/2" (9cm)

NO IMAGE
AVAILABLE

P29181 : TUPPER TYPE BALL CHAINS W/HOOKS

Pkg. of 4, 21cm long

NO IMAGE
AVAILABLE

KM48350 : SCREWDRIVER

2.5mm hexagonal dia., 9 7/8" (25cm)

NO IMAGE
AVAILABLE

KM52782 : JOHNSON RETRACTOR

5" (12.5cm)

NO IMAGE
AVAILABLE

KM48356 : SCREWDRIVER

6mm dia., 10" (25.5cm), Cruciform

NO IMAGE
AVAILABLE

KM52122 : Meyering Finger Retractor #3

4mm width x 8mm depth blade, with teeth, 7" (18.0cm)

NO IMAGE
AVAILABLE

KM41164 : HEGAR BAUMGARTNER NEEDLE HOLDER

Serrated, 5 3/4" (14.5cm)

NO IMAGE
AVAILABLE

KM47950B : SCREWDRIVER MINI HEX

1.5mm x 2mm tip, 2" (5cm)

NO IMAGE
AVAILABLE

KM41068 : WEBSTER NEEDLE HOLDER

Smooth, 5" (12.5cm)

NO IMAGE
AVAILABLE

KM41300 : MAYO HEGAR NEEDLE HOLDER

5" (12.5cm), TC

NO IMAGE
AVAILABLE

KM35380 : Super-Cut Metzenbaum Scissors

Straight, 6" (15.0cm)

NO IMAGE
AVAILABLE

KM51164 : Meyering Laminectomy Retractor, Large

sharp teeth, 50mm width blade, 10" (25.5cm)

NO IMAGE
AVAILABLE

P23587 : Browne Deltoid Retractor

60 x 65mm blade, 10 3/4" (27.5cm)

NO IMAGE
AVAILABLE

380200 : Davis Brain Spatula

380200: Davis Brain Spatula, 1/4" (6.0mm) wide blade, 7-1/8" (18cm) overall length.

NO IMAGE
AVAILABLE

KM47134 : BONE REDUCTION FORCEPS

Pointed tips, 8" (21 cm)

NO IMAGE
AVAILABLE

KM52266 : Volkman Rake Retractor

Sharp, 6-pronged, 1 1/8" (3.0cm) width x 1/2" (1.5cm) depth blade, 8 1/2" (21.5cm)

NO IMAGE
AVAILABLE

P29025 : BUNNELL DISSECTING PROBE

Angled tip, 5" (12.5cm)

NO IMAGE
AVAILABLE

KM51214 : Richardson-Eastmann Retractor, Large

10 5/8" (27.0cm)

NO IMAGE
AVAILABLE

KM52428 : VICKERS UNIVERSAL RETRACTOR

With two blunt prongs to secure placement. - 2 x 2 blunt prongs, - 8 mm deep x 17 mm high, - maximum opening 30 mm, - center blade 10 mm x 16 mm, - 2-1/2" x 3" (64 x 75 mm)

AVAILABLE

NO IMAGE
AVAILABLE

KM48214 : Aufranc Cobra Retractor

hollow handle, smooth pointed 10mm tip, 1 1/4" blade width, 10 1/2" (26.5cm)

NO IMAGE
AVAILABLE

KM52769 : GELPI RETRACTOR

with ball stops, 7" (18cm)

NO IMAGE
AVAILABLE

KM55410 : VERTEBRAL BODY DISTRACTOR, SMALL

NO IMAGE
AVAILABLE

KM51755 : Ribbon Retractor

1" (25mm) width blade, 7 1/2" (19.0cm)

NO IMAGE
AVAILABLE

P37417 : Deltoid Retractor

40mm x 50mm, 10 3/4" (27.5cm)