

## Plastic

NO IMAGE  
AVAILABLE

**065202 : KLEINERT-KUTZ SYNOVECTOMY RONGEUR**

Single action, Curved, 5 1/4" (13.5cm)

NO IMAGE  
AVAILABLE

**065212 : KLEINERT-KUTZ RONGEUR**

Double action, 10 lightly curved, cup 18.5 x 5mm, 7 1/4" (18.5cm)

NO IMAGE  
AVAILABLE

**065300 : KLEINERT-KUTZ ELEVATOR AND DISSECTOR**

Double-ended, Blade A 2.4mm Blade B 3.2mm, 7 3/4" (19.5cm)

NO IMAGE  
AVAILABLE

**KM44224 : MICRO FORCEPS**

Curved, tying platform , 1.2 mm tip, 6" (15 cm)

NO IMAGE  
AVAILABLE

**KM44226 : MICRO FORCEPS**

Curved, tying platform , 1.2 mm tip, 7" (18 cm)

NO IMAGE  
AVAILABLE

**KM44227 : MICRO FORCEPS**

Straight, tying platform, 1.2 mm tip, 7" (18 cm)

NO IMAGE  
AVAILABLE

**KM69459 : EXCAVATOR W/O HOLE # 57-0 1.0MM CURETTE**

5" (12.5 cm)

NO IMAGE  
AVAILABLE

**KM69461 : EXCAVATOR W/O HOLE # 57-1 1.5MM CURETTE**

5" (12.5 cm)

NO IMAGE  
AVAILABLE

**KM69465 : EXCAVATOR W/O HOLE # 57-3 2.5MM CURETTE**

5" (12.5 cm)

NO IMAGE  
AVAILABLE

**KM69478 : EXCAVATOR W/ HOLE # 58-1 1.5 MM CURETTE**

5" (12.5 cm)

NO IMAGE  
AVAILABLE

**KM69480 : EXCAVATOR W/ HOLE # 58-2 2.0 MM CURETTE**

5" (12.5 cm)

NO IMAGE  
AVAILABLE

**KM69484 : EXCAVATOR W/ HOLE # 58-3 2.5 MM CURETTE**

5" (12.5 cm)

NO IMAGE  
AVAILABLE

**P18591 : MICRO VAS CLAMP**

Single 0.6 - 1.5 mm



**P18593 : SINGLE CLAMP**

1.0 mm - 2.25 mm



**P18604 : Approximator Clamp**

0.6 - 1.5mm, S



**P18606 : Approximator Clamp**

1.0 - 2.25mm, A, w/o frame



**P18629 : Artery Approximator Clamp**

Double, 0.6 - 1.5mm, B



**P18630 : ARTERY APPROXIMATOR CLAMP**

Double, 0.6 - 1.5mm, B



**P18633 : Artery Approximator Clamp**

Double, 1.0 - 2.25mm, B



**P18655 : Hand Applied Single Clamps**

One pair

AVAILABLE

NO IMAGE  
AVAILABLE

**P18659 : Approximator Clamp**

W/o frame

NO IMAGE  
AVAILABLE

**P23454 : BONE RASP #64**

Double-ended, 3 mm / 4 mm, 7" (18 cm)

NO IMAGE  
AVAILABLE

**P25642 : JETER-VAN SICKLE SAGITTAL BONE CLAMP**

Heavy, 7" (18 cm)

NO IMAGE  
AVAILABLE

**P25655 : Dunn Dautrey Type Anterior Condylar Retractor**

6.5 x 60mm

NO IMAGE  
AVAILABLE

**P25657 : Dunn Dautrey Type Super Condylar Retractor**

6.5 x 50mm