

Home // Cholangiography // Cholangiography Catheters // TAUT® CHOLANGIOGRAM CATHETER

Catalog # :

~~18400~~

TAUT® CHOLANGIOGRAM CATHETER

Description:

Catheter 4.5 Fr. / 1.5 mm dia. Tip 4.5 Fr. / 1.5 mm dia. Length 30 in. / 76 cm. With metal support tube., 10/BX

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699