

Home // Cholangiography // Cholangiography Catheters // TAUT® CHOLANGIOGRAM CATHETER

Catalog # : 20018-

M57

## TAUT® CHOLANGIOGRAM CATHETER

### Description:

Catheter 4.5 Fr. / 1.5 mm dia. Tip 6.5 Fr. / 2.2 mm dia. Length 18 in. / 45.7 cm. With metal support tube., 10/BX

NO IMAGE  
AVAILABLE

United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**