

Home // Ligation Systems // Weck® Hem-o-lok® Polymer Ligation System - Appliers // Weck® Hem-o-lok® Endo5® Ligating Clip Appliers

Catalog # :

~~544965L~~

Weck® Hem-o-lok® Endo5® Ligating Clip Appliers

Description:

Medium-Large, 17 3/4" (45 cm) applier, 5 mm

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699